

Indicating Desiccant Packets

These Indicating Desiccant Packets provide reliable moisture protection with the added benefit of visual monitoring. These packets contain color-changing indicator beads that signal when the desiccant has reached its moisture absorption capacity, making it easy to know when replacement is needed.



FS404

SPECIFICATIONS

Construction

- Each packets contains (1) gram of cobaltous silica gel with millions of pores
- Adheres to FDA guidelines

Quantity

- 1,000 per case

Dimensions

- 1.75" L X 0.75" W X 0.25" H

Performance

- Soaks up water vapor
- Eliminates moisture and dampness
- Gel turns pink when moisture is absorbed

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI