

Fingernail Clippings Evidence Collection Kit

Fingernail Clippings Evidence Collection Kit is designed for the safe and effective collection of nail evidence in assault investigations. Used when fingernails may break during a struggle, this kit helps recover clippings for comparison and potential DNA analysis.



FS460

SPECIFICATIONS

Construction

- (1) Tamper-evident, pre-sealed kit envelope
- (1) Kit instruction sheet
- (1) Metal fingernail clipper
- (1) Right-hand envelope containing one paper bundle
- (1) Left-hand envelope containing one paper bundle
- (1) Tamper-evident police evidence seal for resealing kit envelope after specimen collection

Dimensions

- Envelope: 5.5" L X 7.5" W

Quantity

- 25 per case

Performance

- Everything needed for fingernail collection for comparing with evidence

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI