


# Fingernail and Fingertip Swabbing Evidence Collection Kit

This Fingernail and Fingertip Swabbing Evidence Collection Kit is designed for cases where a victim reports scratching an assailant but no visible debris is present under the nails. Using sterile, moistened swabs, the kit allows investigators to collect trace DNA from beneath fingernails and fingertips for laboratory analysis and comparison.



FS462

## SPECIFICATIONS

### Construction

- (1) Tamper-evident, pre-sealed Kit Envelope
- (1) Kit Instruction Sheet
- (2) Sterile, cotton-tipped Swabs (2/pkg)
- (1) 5 cc Sterile Water
- (1) Right-Hand Swab Box
- (1) Left-Hand Swab Box
- (1) Tamper-evident Police Evidence Seal for resealing kit envelope after specimen collection

### Performance

- For collection of trace DNA from under fingernails

### Dimensions

- Envelope: 5.5" X 7.5"

### Quantity

- 25 per case

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI