

Glassine Weight and Trace Paper

Glassine Weight and Trace Paper is designed for the secure handling and preservation of delicate trace evidence. Made from smooth, semi-transparent glassine, this paper resists moisture, air, and static, making it ideal for collecting and transferring materials like hairs, fibers, paint chips, and other micro evidence without loss or contamination.



SPECIFICATIONS

Construction

- Provided in a re-closable box
- Smooth, lightweight, and dust-free
- Strong enough for analytical samples

Dimensions

- Small: 3" L X 3" W
- Medium: 4" L X 4" W
- Large: 6" L X 6" W

Quantity

- 500 sheets per box

Variations

Model No.	Description
FS466	Small Glassine Weight and Trace Paper
FS467	Medium Glassine Weight and Trace Paper
FS468	Large Glassine Weight and Trace Paper

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI