

P-DMAC Paper

P-DMAC Paper is designed for developing latent fingerprints on porous, heat-sensitive surfaces such as thermal paper and receipts. It reacts with fingerprint residues to produce fluorescent ridge detail when viewed under alternate light sources.



FS470

SPECIFICATIONS

Performance

- 6 month shelf life
- Best results usually visualized between 475nm and 515nm with an orange barrier filter or glasses
- To extend longevity of this product, always store in the provided amber sleeve.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI