


GSR Collection Kit, 12 Swabs

This GSR Collection Kit is designed for the efficient collection of gunshot residue from a subject's hands and other surfaces. Each kit includes 12 sterile swabs for thorough sampling, helping ensure sufficient material is collected for laboratory analysis. Packaged with tamper-evident components, it supports proper evidence handling, preservation, and chain of custody.



FS479

SPECIFICATIONS

Construction

- (1) Tamper-evident, pre-sealed kit envelope
- Kit instruction sheet
- (1) Gunshot residue analysis information form
- (1) Pair of nitrile barrier gloves
- (1) Resealable bag containing (12) cotton-tipped, plastic shaft swabs
- (1) 3 mL dispenser containing a 5% Nitric acid Solution
- (1) Control tube and cap
- (1) Cartridge case tube and cap
- (1) Right back, right palm, left back, and left palm tubes and caps
- (1) Tamper-evident police evidence seal for resealing kit after evidence collection

Dimensions

- Envelope: 6" L X 9" W
- Resealable Bag: 2" L X 3" W

Quantity

- 10 per case

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI