

# GSR Collection Kit, 2 Stubs

These kits are designed for agencies using AA analysis to collect gunshot residue samples. It includes swabs, nitric acid solution, labeled collection tubes, and all necessary supplies for secure, tamper-evident evidence collection.



FS480

## SPECIFICATIONS

### Construction

- (1) Tamper-evident, pre-sealed Kit Envelope
- Kit instruction sheet
- (1) Gunshot residue analysis information form
- (1) Pair of disposable plastic gloves
- (1) Right hand vial with cap containing an SEM stub covered with double-sided adhesive carbon tape
- (1) Left hand vial with cap containing an SEM stub covered with double-sided adhesive carbon tape
- (1) Tamper-evident police evidence seal for resealing kit after evidence collection

### Dimensions

- Envelope: 6" L X 9" W

### Quantity

- 10 per case

### Performance

- For agencies whose crime labs perform GSR analysis by SEM on two collection stubs

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI