


GSR Collection Kit, 4 Stubs, 12 Swabs

GSR Collection Kit, 4 Stubs, 12 Swabs is designed for the thorough collection of gunshot residue from a subject's hands and nearby surfaces. The kit includes four adhesive stubs for capturing particulate evidence and twelve sterile swabs for additional sampling, providing a complete approach to GSR collection.



FS481

SPECIFICATIONS

Construction

- (1) Tamper-evident, pre-sealed kit envelop
- Kit instruction sheet
- (1) Gunshot residue analysis information form
- (1) Pair of disposable plastic gloves
- (1) Resealable bag containing
- (12) Cotton-tipped, plastic shaft swabs
- (1) 3 mL dispenser containing a 5% nitric acid solution
- (1) Control tube and cap
- (1) Cartridge case tube and cap
- (1) Right back, right palm, left back, and left palm tubes and caps
- (1) Right back, right palm, left back, and left palm vials with caps containing SEM stubs covered with double-sided adhesive carbon tape
- (1) Tamper-evident police evidence seal for resealing kit after evidence collection

Quantity

- 10 per case

Dimensions

- Envelope: 6" L X 9" W
- Resealable bag: 2" L X 3" W

Performance

- Perform quick screen using AA and then more definitive analysis by SEM

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI