

GSR Collection Kit, 4 Stubs

These kits are designed for SEM-based gunshot residue collection, capturing Right Back, Right Palm, Left Back, and Left Palm hand samples. It includes tamper-evident, sealed components with adhesive carbon tape stubs for secure evidence collection.



FS482

SPECIFICATIONS

Construction

- (1) Tamper-evident, pre-sealed kit envelope
- Kit instruction sheet
- (1) Gunshot residue analysis information form
- (1) Pair of disposable plastic gloves
- (1) Right and left back vials with cap containing a SEM stub covered with double-sided adhesive carbon tape
- (1) Right and left palm vials with cap containing a SEM stub covered with double-sided adhesive carbon tape
- (1) Tamper-evident police evidence seal for resealing kit after evidence collection

Dimensions

- Envelope: 6" L X 9" W

Quantity

- 10 per case

Performance

- For law enforcement agencies whose crime labs use SEM as their method of GSR analysis request "Right Back", "Left Back", "Right Palm", and "Left Palm" samples be collected from the suspect's hands

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI