

Surface GSR Collection Kit

These kits are designed for collecting gunshot residue from surfaces such as tabletops, vehicles, doorframes, or clothing. It includes labeled SEM collection vials with adhesive carbon stubs and all necessary materials for clear, traceable evidence collection.



FS483

SPECIFICATIONS

Construction

- (1) Tamper-evident, pre-sealed kit envelope
- Kit instruction sheet
- (1) Surface gunshot residue analysis information form
- (1) Pair of disposable plastic gloves
- (1) Sample numbers 1-4 vials with caps containing SEM stubs covered with double-sided adhesive carbon tape
- (1) Tamper-evident police evidence seal for resealing kit after evidence collection

Performance

- Collects GSR from surfaces other than the suspect's hand

Dimensions

- Envelope: 6" L X 9" W

Quantity

- 10 per case

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI