

SEM Collection Devices

SEM collection devices are used for collecting gunshot residue and other trace evidence for SEM analysis. Each vial contains an adhesive carbon stub with a write-on label for sample identification. Available in Standard, Q-Stub, or D-Stub formats.



SPECIFICATIONS

Construction

- Each vial contains a SEM stub covered with a double-sided adhesive carbon tape
- Write-on label

Performance

- For collection of explosive residue collection, inorganic particle matter, preparing standards for future comparison, paint chips, glass, soil, fiber, or hair

Materials

- Plastic
- Carbon tape

Quantity

- 50 per case

Variations

Model No.	Description
FS484	Standard SEM Collection
FS485	D-Stub SEM Collection
FS486	Q-Stub SEM Collection

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI