



Straight Stainless Steel Hemostat

Surgical-grade stainless steel hemostats are used for handling evidence, holding documents, and other general lab tasks. They are durable, easy to clean with alcohol or other cleaners, and can be autoclaved for sterilization. Available in three different sizes.



SPECIFICATIONS

Performance

- For picking up evidence or holding documents
- Surgical grade
- Easily cleaned with alcohol or other cleaners
- Can be autoclaved

Dimensions

- FS487: 5.5" L
- FS488: 8" L
- FS489: 12" L
- FS490: 16" L

Materials

- Stainless steel

Variations

Model No.	Description
FS487	Straight Hemostat, 5.5"
FS488	Straight Hemostat, 8"
FS489	Straight Hemostat, 12"
FS490	Straight Hemostat, 16"

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI