

Nickel Plated Forceps

Nickel Plated Forceps are designed for precise handling of small and delicate evidence items in forensic and laboratory settings. Constructed from durable nickel-plated steel, they offer strength, corrosion resistance, and a secure grip for controlled manipulation.




FS495



FS497



FS496



FS498

SPECIFICATIONS

Performance

- Easily cleaned with alcohol and other solvents
- Autoclavable
- Ideal for picking up trace evidence such as hair, glass shards, etc

Materials

- Nickel

Variations

Model No.	Description
FS495	Straight Forceps
FS496	Dumont Forceps
FS497	Carmalt Forceps
FS498	Angle Point Forceps

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI