


16 Gauge Needle

The 16 Gauge Needles are designed for precise puncturing and sample handling in forensic and laboratory applications. Constructed from durable stainless steel, it provides strength and consistency for controlled use. 100 needles per pack.



FS499

SPECIFICATIONS

Construction

- Sterile
- 16 gauge
- Hypodermic needles

Dimensions

- 4" Length

Materials

- Stainless steel

Quantity

- 100 per pack

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI