

# Aluminum Fuming Trays

Aluminum Fuming Trays are designed for use in cyanoacrylate fuming processes to develop latent fingerprints. Constructed from lightweight, heat-conductive aluminum, these trays provide an even surface for controlled heating and vapor distribution.



FS501

## SPECIFICATIONS

### Performance

- Accelerate the cyanoacrylate fuming process by using a hot plate (sold separately) and aluminum dish inside a fuming chamber

### Quantity

- 100 per case

### Option and Accessories

- FS502: Hot Plate, 120v

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI