

# Peel & Fume Cyanoacrylate Latent Print Fuming Pouches

Peel & Fume pouches are a one-step cyanoacrylate system for rapid latent print development in a fuming chamber. Each pouch contains gel superglue that activates on exposure, developing prints in 10–20 minutes.



## SPECIFICATIONS

### Construction

- Peel-open pouch
- Easy to use, one-part fuming system

### Materials

- Foil

### Performance

- Develops in 10–20 minutes

### Quantity

- FS510: 20 per package
- FS511: 100 per package

## Variations

Model No.	Quantity
FS510	Fuming Pouches, 20 per package
FS511	Fuming Pouches, 100 per package

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI