

Hemident Presumptive Blood ID Kit

Hemident is a nondestructive presumptive test for bloodstains that preserves evidence for further laboratory analysis. It provides rapid field results, producing a blue-green color change in the presence of mammal blood.



FS512

SPECIFICATIONS

Performance

- For presumptive identification of bloodstain
- Non-destructive and enables the investigator to retain the evidence intact for further analysis
- Reagents are capable of identifying down to one part per million of blood

Construction

- (2) Ampoules within a plastic pouch

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI