

Presumptive Blood ID Kit

This Presumptive Blood ID Kit is a rapid presumptive blood test that distinguishes human blood from animal blood in 2–3 minutes. A simple field procedure produces a visual result with two blue bars confirming human blood, helping prioritize evidence collection and reduce unnecessary lab testing.



SPECIFICATIONS

Construction

- Hexagon OBTI reagent bottle
- Sample collection stick

Quantity

- FS513: 6 per box
- FS514: 24 per box

Performance

- Distinguishes between human and animal blood
- Results in 2–3 minutes
- If the sample is human blood, two blue bars will appear

Variations

Model No.	Description
FS513	Presumptive Blood ID Kit, 6 per box
FS514	Presumptive Blood ID Kit, 24 per box

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI