

Fingerprint Rehydration Solution

This Fingerprint Rehydration Solution improves fingerprint ridge detail for both living and deceased individuals using a non-toxic spray application. It enhances print clarity by rehydrating skin tissue for faster, more defined fingerprint impressions.



SPECIFICATIONS

Construction

- Clear, odorless, non-toxic, and completely non-volatile solution
- Injectable tissue rehydration fluid

Performance

- Rub the solution over the intended ridge detail and sufficient time should be allowed for the solution to absorb into the skin tissue

Capacity

- FS515: 0.25 oz
- FS516: 1 oz
- FS517: 2 oz
- FS518: 4 oz

Variations

| Model No. | Description |
|-----------|-------------------------------|
| FS515 | Rehydration Solution, 0.25 oz |
| FS516 | Rehydration Solution, 1 oz |
| FS517 | Rehydration Solution, 2 oz |
| FS518 | Rehydration Solution, 4 oz |

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI