

Ninhydrin Spray, Various Sizes

This quick-drying Ninhydrin Spray is an easy-to-use, extremely sensitive method for developing latent prints on paper and other surfaces. It reacts with proteins and amino acids in fingerprints to produce a purple print when heated.



SPECIFICATIONS

Construction

- Pump misting bottle
- Detailed instruction sheet is supplied

Performance

- Reacts to the proteins and amino acids found in fingerprints

Options & Accessories

- Regular or extra-strength solutions available
- 8 or 16 oz options available

Variations

Model No.	Description
FS522	Regular Ninhydrin Spray, 8 oz
FS523	Regular Ninhydrin Spray, 16 oz
FS524	Extra-Strength Ninhydrin Spray, 8 oz
FS525	Extra-Strength Ninhydrin Spray, 16 oz

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI