

Ninhydrin Forensic Solution

Ninhydrin is an easy and effective method for developing latent prints on paper and is widely used by latent print examiners. It can be mixed with 3M Novec Engineered Fluid HFE-7100 for a less toxic and non-flammable solution.



SPECIFICATIONS

Construction

- Working solution has a one-year shelf life

Materials

- Ninhydrin
- Ninhydrin ethanol formula
- HFE-7100 formula

Performance

- Easily develop latent print on paper
- HFE-7100 formula can be mixed with the Ninhydrin ethanol formula for a less toxic and non-flammable solution

Variations

Model No.	Description	Capacity
FS528	Ninhydrin Powder	25 grams
FS529	Ninhydrin Ethanol Formula	250 ml
FS530	Ninhydrin Ethanol Formula	500 ml
FS531	Ninhydrin HFE-7100 Formula	250 ml

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI