

# Ninhydrin Forensic Solution

Ninhydrin is an easy and effective method for developing latent prints on paper and is widely used by latent print examiners. It can be mixed with 3M Novec Engineered Fluid HFE-7100 for a less toxic and non-flammable solution.



## SPECIFICATIONS

### Construction

- Crystals have an indefinite shelf life
- Working solution had one year shelf life

### Performance

- Easily develop latent print on paper
- HFE-7100 formula can be mixed with ninhydrin ethanol formula for a less toxic and non-flammable solution

## Variations

Model No.	Description
FS528	Ninhydrin Powder, 25 grams
FS529	Ninhydrin Ethanol Formula, 250 mL
FS530	Ninhydrin Ethanol Formula, 500 mL
FS531	Ninhydrin HFE-7100 Formula, 250 mL

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI