

Fingerstick Lancets

Fingerstick Lancets are designed for quick, sterile, and minimally invasive blood sample collection. Engineered for single-use safety, each lancet delivers a controlled puncture for consistent capillary blood flow, making them ideal for forensic, medical, and laboratory applications.



FS540

SPECIFICATIONS

Construction

- Protective lancet base
- 21 gauge stick

Performance

- Twist off base, hold against subject's finger, and press the green button

Quantity

- 100 per box

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI