

Activator Solution

The Activator Solution is used to activate the Fluorescent Blood Lifting Strips prior to lifting bloody impression evidence for analysis. It is supplied in an 8 oz bottle and is used as part of a quick and easy process for enhancing forensic impressions.



FS544

SPECIFICATIONS

Performance

- Activates the fluorescent blood lifting strip (sold separately)
- Provides a safe, non-toxic fluorogenic method for lifting bloody impression evidence

Capacity

- 8 oz

Options and Accessories

- FS612: 2" X 2" Fluorescent Blood Lifting Strips, (8) per package
- FS613: 2" X 4" Fluorescent Blood Lifting Strips, (4) per package
- FS614: 4" X 4" Fluorescent Blood Lifting Strips, (2) per package
- FS615: 8" X 8" Fluorescent Blood Lifting Strips, (1) per package

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI