

Luminol Blood Detection Reagent Spray

Luminol Blood Detection Reagent Spray is used to detect trace or invisible blood at crime scenes. When applied in a darkened room, it produces a white/blue glow where blood is present, allowing the pattern to be photographed.



SPECIFICATIONS

Construction

- Nitrile barrier gloves
- Vial of activation fluid

Materials

- Nitrile
- Luminol

Performance

- Luminol glows blue/white when it comes into contact with blood in a darkened room

Capacity

- FS545: 4 oz
- FS546: 8 oz
- FS547: 16 oz

Variations

Model No.	Description
FS545	Luminol Blood Detection Reagent Spray, 4 oz
FS546	Luminol Blood Detection Reagent Spray, 8 oz
FS547	Luminol Blood Detection Reagent Spray, 16 oz

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI