

Hydro-Print Wet Latent Print Developer

The Hydro-Print Wet Latent Print Developer is formulated for developing latent prints on wet or damp surfaces such as glass, plastic, and metal. It adheres to fingerprint oils for easy visualization and can be photographed and lifted after development, with a simple mix-and-spray application.



SPECIFICATIONS

Construction

- Hydro-print (small particle reagent) bottle comes with:
- 1 detailed instruction sheet
- Specifically formulated for developing latent prints found outdoors on wet or damp surfaces
- Most effective if used within 6 months of manufactured date
- Store in a cool, dark place

Performance

- Add about 16 oz of water to the spray bottle, then shake the bottle until it is homogeneous
- Adheres to the oils found in fingerprints

Capacity

- 16 oz

Variations

Model No.	Description	Color
FS556	Latent Print Developer	Black
FS557	Latent Print Developer	White

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI