

Hydro-Print Wet Latent Print Developer

Hydro-Print (Small Particle Reagent) is formulated for developing latent prints on wet or damp surfaces such as glass, plastic, and metal. It adheres to fingerprint oils for easy visualization and can be photographed and lifted after development, with a simple mix-and-spray application.



SPECIFICATIONS

Construction

- Bottle contains small particle reagent
- Stored in a cool, dark place

Performance

- For developing prints found outdoors on wet or damp surfaces
- Adheres to the oils found in fingerprints

Capacity

- 16 oz

Options & Accessories

- Black or White option available

Variations

Model No.	Description
FS556	Black Latent Print Developer
FS557	White Latent Print Developer

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI