

# Trace Metal Detection Spray

Trace Metal Detection Spray is a forensic reagent designed to reveal trace metal residues left on a suspect's hands after contact with a weapon or metal object. Applied directly to the skin and viewed under short wave ultraviolet light, the spray helps investigators identify telltale patterns and impressions of metallic items.



FS563

## SPECIFICATIONS

### Construction

- Spray bottle filled with reagent

### Performance

- When sprayed on the suspect's hand, short wave UV light reveals traces of a metal tool or weapon

### Capacity

- 8 oz bottle

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI