

# Deluxe Hinged Print Lifter Kit

The Deluxe Hinged Print Lifter Kit contains 144 hinged lifters in multiple sizes and backing colors for preserving developed latent prints. It uses crystal clear adhesive mylar and vinyl backing cards to securely lift, document, and store fingerprint evidence.



FS566

## SPECIFICATIONS

### Construction

- Hinged Print Lifters: White, black, and transparent in 1.5" L X 2" W, 2" L X 4" W, and 4" L X 4" W (16 of each combination)

### Performance

- Pertinent information pertaining to the print can be recorded on the reverse side of the lifter
- Peel off the liner to expose adhesive, lift the print, then close the Mylar to permanently secure it

### Dimensions

- Box: 7" L X 5" W X 1.875" H
- Covers: 0.005" thick
- Back Cards: 0.01" thick

### Materials

- Mylar
- Vinyl

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI