

Basic Hinged Print Lifter Kit

This Hinged Print Lifter Kit is a compact plastic case containing 48 lifters with white, black, or transparent backing cards, designed for efficiently collecting developed fingerprints. Each lifter uses a clear adhesive mylar cover hinged to a vinyl backing card, allowing users to peel, lift, and securely preserve prints with ease.



SPECIFICATIONS

Construction

- (24) 1.5" X 2" Hinged Print Lifters
- (12) 2" X 4" Hinged Print Lifters
- (12) 4" X 4" Hinged Print Lifters

Quantity

- 48 per box

Materials

- Mylar
- Vinyl

Performance

- Lifts developed latent fingerprints from surfaces
- Preserves prints for analysis and evidence
- Provides clear, secure transfer using adhesive lifters
- Supports long-term storage and examination

Dimensions

- Box: 4.75" L X 4.75" W X 1.375" H
- Covers: 0.005" thick
- Back Cards: 0.01" thick

Options & Accessories

- Available in three different colors: White, black, or transparent

Variations

Model No.	Description
FS567	Black Hinged Print Lifter Kit
FS568	White Hinged Print Lifter Kit
FS569	Transparent Hinged Print Lifter Kit

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI