

Hinged Print Lifters - 1.5" X 2"

Crystal-clear hinged latent print lifters designed for secure collection, preservation, and documentation of developed fingerprints. Constructed with adhesive-backed Mylar and durable vinyl backing cards to help protect and preserve lifted prints during transport and examination.



SPECIFICATIONS

Construction

- Clear adhesive backing
- Release liner

Quantity

- 24 per box

Materials

- Mylar
- Vinyl

Performance

- Lifts developed latent prints from surfaces
- Seals prints for permanent preservation
- Allows written case notes on backing card
- Used for fingerprints, palm prints, and footprint residue
- Provides clear visual preservation of evidence

Dimensions

- 1.5" X 2"
- .005" thick adhesive-backed mylar
- .01" thick vinyl backing cards

Variations

Model No.	Description	Color
FS570	Hinged Print Lifters	Black
FS571	Hinged Print Lifters	White
FS572	Hinged Print Lifters	Transparent

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI