

Hinged Print Lifters - 2" X 4"

Crystal-clear hinged latent print lifters designed for secure collection, preservation, and documentation of developed fingerprints. Constructed with adhesive-backed Mylar and durable vinyl backing cards to help protect and preserve lifted prints during transport and examination.



SPECIFICATIONS

Construction

- Clear adhesive backing
- Release liner

Quantity

- 12 per box

Materials

- Mylar
- Vinyl

Performance

- Lifts developed latent prints from surfaces
- Seals prints for permanent preservation
- Allows written case notes on backing card
- Used for fingerprints, palm prints, and footprint residue
- Provides clear visual preservation of evidence

Dimensions

- 2" X 4"
- .005" thick adhesive-backed mylar
- .01" thick vinyl backing cards

Variations

| Model No. | Description | Color |
|-----------|----------------------|-------------|
| FS573 | Hinged Print Lifters | Black |
| FS574 | Hinged Print Lifters | White |
| FS575 | Hinged Print Lifters | Transparent |

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI