

Deluxe Hinged Print Lifter Kit, Scaled

The Deluxe Scaled Hinged Print Lifter Kit includes 144 lifters in multiple sizes and backing colors, stored in a durable plastic case. Each lifter features clear adhesive mylar hinged to a vinyl backing card with printed scales for accurate fingerprint collection and documentation.



FS579

SPECIFICATIONS

Construction

- Scaled Hinged Print Lifters: White, black, and transparent in 1.5" L X 2" W, 2" L X 4" W, and 4" L X 4" W (16 of each combination)

Quantity

- 144 per box

Performance

- Lifts and preserves developed fingerprints with measurement reference
- Provides clear, secure transfer using adhesive mylar
- Supports accurate documentation with built-in scale
- Allows labeling on reverse side for evidence tracking

Dimensions

- Box: 7" L X 5" W X 1.875" H
- Covers: 0.005" thick
- Back Cards: 0.01" thick

Materials

- Mylar
- Vinyl

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI