

Basic Scaled Hinged Print Lifter Kit

This Hinged Print Lifter Kit is a compact plastic case containing 48 lifters with white, black, or transparent backing cards, designed for efficiently collecting developed fingerprints. Each lifter uses a clear adhesive mylar cover hinged to a vinyl backing card, allowing users to peel, lift, and securely preserve prints with ease.



SPECIFICATIONS

Construction

- (24) 1.5" X 2" Scaled Hinged Print Lifters
- (12) 2" X 4" Scaled Hinged Print Lifters
- (12) 4" X 4" Scaled Hinged Print Lifters

Quantity

- 48 per box

Performance

- Lifts latent fingerprints from surfaces
- Preserves prints on clear film for analysis
- Secures prints onto backing cards for storage
- Includes scale for size reference
- Allows notes for documentation

Dimensions

- Box: 4.75" L X 4.75" W X 1.375" H
- Covers: 0.005" thick
- Back Cards: 0.01" thick
- Scale: measures 40 mm X 23 mm

Materials

- Mylar
- Vinyl

Variations

Model No.	Description
FS580	Black Scaled Hinged Print Lifter Kit
FS581	White Scaled Hinged Print Lifter Kit
FS582	Transparent Scaled Hinged Print Lifter Kit

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI