

# Scaled Hinged Print Lifters - 1.5" X 2"

A premium, hinged adhesive lifter engineered specifically to lift, collect, and preserve latent fingerprints from a variety of surfaces. Each lifter features a built-in, precise measurement scale for immediate forensic documentation, along with a durable vinyl backing card designed for long-term, secure evidence storage and case labeling.



## SPECIFICATIONS

### Construction

- Clear adhesive backing
- Release liner
- 40 mm x 23 mm scale on bottom and right side

### Quantity

- 24 per box

### Performance

- Lifts developed latent fingerprints from surfaces
- Preserves prints on backing card for evidence storage
- Provides built-in scale for measurement reference
- Allows case notes on reverse side of backing card

### Dimensions

- 1.5" X 2"
- .005" thick adhesive-backed mylar
- .01" thick vinyl backing cards

### Materials

- Mylar
- Vinyl

## Variations

Model No.	Description	Color
FS583	Scaled Hinged Print Lifters	Black
FS584	Scaled Hinged Print Lifters	White
FS585	Scaled Hinged Print Lifters	Transparent

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI