

Scaled Hinged Print Lifters - 2" X 4"

A premium, hinged adhesive lifter engineered specifically to lift, collect, and preserve latent fingerprints from a variety of surfaces. Each lifter features a built-in, precise measurement scale for immediate forensic documentation, along with a durable vinyl backing card designed for long-term, secure evidence storage and case labeling.



SPECIFICATIONS

Construction

- Clear adhesive backing
- Release liner
- 40 mm x 86 mm scale on bottom and right side

Quantity

- 12 per box

Performance

- Lifts developed latent fingerprints from surfaces
- Preserves prints on backing card for evidence storage
- Provides built-in scale for measurement reference
- Allows case notes on reverse side of backing card

Dimensions

- 2" X 4"
- .005" thick adhesive-backed mylar
- .01" thick vinyl backing cards

Materials

- Mylar
- Vinyl

Variations

Model No.	Description	Color
FS586	Scaled Hinged Print Lifters	Black
FS587	Scaled Hinged Print Lifters	White
FS588	Scaled Hinged Print Lifters	Transparent

MOPEC
GROUP

mopec
ELEVATING ANTHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI