

Scaled Hinged Print Lifters - 4" X 4"

A hinged adhesive lifter used to collect and preserve latent fingerprints with an integrated measurement scale and vinyl backing card for documentation and storage.



SPECIFICATIONS

Construction

- Clear adhesive backing
- Release liner
- 90 mm x 80 mm scale on bottom and right side

Quantity

- 24 per box

Performance

- Lifts developed latent fingerprints from surfaces
- Preserves prints on backing card for evidence storage
- Provides built-in scale for measurement reference
- Allows case notes on reverse side of backing card

Dimensions

- 4" X 4"
- .005" thick adhesive-backed mylar
- .01" thick vinyl backing cards

Materials

- Mylar
- Vinyl

Variations

Model No.	Description	Color
FS589	Scaled Hinged Print Lifters	Black
FS590	Scaled Hinged Print Lifters	White
FS591	Scaled Hinged Print Lifters	Transparent

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI