

# Hinged Palm Print Lifter

Hinged Palm Print Lifters are designed for lifting and preserving developed palm prints using a clear adhesive-backed mylar cover. Once applied, the print is secured onto a vinyl backing card for permanent storage and documentation.



## SPECIFICATIONS

### Construction

- Clear adhesive backing
- Release liner

### Quantity

- 12 per box

### Materials

- Mylar
- Vinyl

### Performance

- Lifts latent fingerprints from surfaces
- Preserves prints on clear film for analysis
- Secures prints onto backing cards for storage
- Includes scale for size reference
- Allows notes for documentation

### Dimensions

- Print: 6" L X 9" W
- Covers: 0.005" thick
- Backing Cards: 0.01" thick

### Options & Accessories

- Available in black, white or transparent

## Variations

Model No.	Description
FS592	Black Hinged Palm Print Lifter
FS593	White Hinged Palm Print Lifter
FS594	Transparent Hinged Palm Print Lifter

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI