

# Atomic Latent Print Powder, Pure White

Atomic Latent Print Powder is a high-performance fingerprint powder formulated to develop clear, sharp latent prints on a variety of surfaces. Its fine particle consistency provides excellent adhesion and contrast, making it ideal for detailed forensic investigations and evidence processing.



\*Photo for reference only. Powder is Pure White.

## SPECIFICATIONS

### Construction

- Smooth, non-caking
- Moisture-free

### Capacity

- Available sizes:  
2 oz - 128 oz

### Performance

- Develops clear latent fingerprints on various surfaces
- Strong adhesion to fingerprint secretions for better detail
- Produces smooth, consistent application without clumping

### Option & Accessories

- More colors available, please reach out to your representative.

## Variations

Model No.	Description	Capacity
FS659	Pure White Latent Print Powder	2 oz
FS660	Pure White Latent Print Powder	4 oz
FS661	Pure White Latent Print Powder	8 oz
FS662	Pure White Latent Print Powder	16 oz
FS663	Pure White Latent Print Powder	32 oz
FS664	Pure White Latent Print Powder	64 oz
FS665	Pure White Latent Print Powder	128 oz

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI