

Magnetic Atomic Latent Print Powders

Magnetic Atomic Latent Print Powders are high-quality magnetic fingerprint powders designed for developing latent prints on challenging nonferrous surfaces. They are used with magnetic applicators to reveal prints without direct brush contact.



SPECIFICATIONS

Construction

- Smooth, non-caking formula and application
- Moisture-free
- Denser weight properties from standard latent print powder

Capacity

- Available in 1 oz, 16 oz or 128 oz jars

Performance

- Develops latent prints on difficult surfaces like glass, wood, leather, and vinyl
- Provides strong contrast for clear visualization of ridge detail
- Works effectively on both smooth and textured materials

Option and Accessories

- Color options:
Regular Black
Midnight Black
White
Gray
- Note: More colors available. Please reach out to your representative.

Variations

Model No.	Description	Capacity
FS694	Regular Black Magnetic Latent Print Powder	1 oz
FS695	Regular Black Magnetic Latent Print Powder	16 oz
FS696	Regular Black Magnetic Latent Print Powder	128 oz
FS697	Midnight Black Magnetic Latent Print Powder	1 oz
FS698	Midnight Black Magnetic Latent Print Powder	16 oz
FS699	Midnight Black Magnetic Latent Print Powder	128 oz
FS700	White Magnetic Latent Print Powder	1 oz
FS701	White Magnetic Latent Print Powder	16 oz
FS702	White Magnetic Latent Print Powder	128 oz
FS703	Gray Magnetic Latent Print Powder	1 oz
FS704	Gray Magnetic Latent Print Powder	16 oz
FS705	Gray Magnetic Latent Print Powder	128 oz

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI