

Magnetic Atomic Latent Print Powders

Magnetic Atomic Latent Print Powders are high-quality magnetic fingerprint powders designed for developing latent prints on challenging nonferrous surfaces. They are used with magnetic applicators to reveal prints without direct brush contact.



SPECIFICATIONS

Construction

- Smooth, non-caking formula and application
- Moisture-free
- Denser weight properties from standard latent print powder

Capacity

- Available in 1 oz, 16 oz or 128 oz jars

Performance

- Develops latent prints on difficult surfaces like glass, wood, leather, and vinyl
- Provides strong contrast for clear visualization of ridge detail
- Works effectively on both smooth and textured materials

Option and Accessories

- Color options:
Red
Silver
Silver/Black
- Note: More colors available. Please reach out to your representative.

Variations

| Model No. | Description | Capacity |
|-----------|---|----------|
| FS706 | Red Magnetic Latent Print Powder | 1 oz |
| FS707 | Red Magnetic Latent Print Powder | 16 oz |
| FS708 | Red Magnetic Latent Print Powder | 128 oz |
| FS709 | Silver Magnetic Latent Print Powder | 1 oz |
| FS710 | Silver Magnetic Latent Print Powder | 16 oz |
| FS711 | Silver Magnetic Latent Print Powder | 128 oz |
| FS712 | Silver/Black Magnetic Latent Print Powder | 1 oz |
| FS713 | Silver/Black Magnetic Latent Print Powder | 16 oz |
| FS714 | Silver/Black Magnetic Latent Print Powder | 128 oz |

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI