

# Dual Purpose Atomic Latent Print Powder

Dual Purpose (Bi-Chromatic) Atomic Latent Print Powders are fingerprint powders designed to adapt to different surface colors for improved visualization of latent prints. They provide strong contrast on both light and dark backgrounds for more effective evidence development.



## SPECIFICATIONS

### Construction

- Smooth, non-caking formula and application
- Moisture-free

### Capacity

- Available in 2 oz, 8 oz or 16 oz jars

### Performance

- Develops clear latent prints on multi-colored surfaces
- Appears light on dark surfaces and dark on light surfaces
- Enhances ridge detail visibility across varied backgrounds

### Options & Accessories

- Available color options:  
Silver/Black  
Silver/Gray  
Silver/Red

## Variations

| Model No. | Description                                   | Capacity |
|-----------|---|----------|
| FS727     | Silver/Black Dual Purpose Latent Print Powder | 1 oz     |
| FS728     | Silver/Black Dual Purpose Latent Print Powder | 1 oz     |
| FS729     | Silver/Black Dual Purpose Latent Print Powder | 1 oz     |
| FS730     | Silver/Gray Dual Purpose Latent Print Powder  | 1 oz     |
| FS731     | Silver/Gray Dual Purpose Latent Print Powder  | 1 oz     |
| FS732     | Silver/Gray Dual Purpose Latent Print Powder  | 1 oz     |
| FS733     | Silver/Red Dual Purpose Latent Print Powder   | 1 oz     |
| FS734     | Silver/Red Dual Purpose Latent Print Powder   | 1 oz     |
| FS735     | Silver/Red Dual Purpose Latent Print Powder   | 1 oz     |

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI