

# Safe-Cracker Atomic Latent Print Powder

Safe-Cracker Latent Print Powder is specially formulated for processing difficult or textured surfaces where traditional powders may be less effective. Designed for enhanced visibility and adhesion, it helps reveal detailed latent prints while minimizing background interference during forensic examinations.



## SPECIFICATIONS

### Construction

- Smooth, non-caking formula and application
- Moisture-free
- Metallic red color

### Capacity

- Available in 2 oz or 16 oz jars

### Performance

- Develops latent prints on metal surfaces
- Enhances contrast on metallic backgrounds
- Produces clear, visible ridge detail
- Effective on items like safes, vaults, and cabinets

## Variations

Model No.	Description	Capacity
FS748	Safe-Cracker Atomic Latent Print Powder	2 oz
FS749	Safe-Cracker Atomic Latent Print Powder	16 oz

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI