


Wet Powder

Wet Powder is a specialized forensic powder designed for developing latent fingerprints on adhesive and non-porous surfaces. Its advanced formulation adheres effectively to moisture-covered evidence, helping investigators recover clear, high-contrast fingerprints in challenging environmental conditions.



FS750



FS751

SPECIFICATIONS

Performance

- Develops fingerprints on adhesive and non-porous surfaces
- Can be used alone or after cyanoacrylate fuming
- Maintains integrity of delicate or sensitive evidence

Capacity

- 250 mL

Option and Accessories

- Available in black (FS750) or white (FS751)

Variations

Model No.	Description
FS750	Black Wet Powder
FS751	White Wet Powder

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI