

Blood Fixative

Blood Fixative is a forensic reagent used to stabilize and preserve fingerprints or impression evidence left in blood before further processing. It helps maintain the clarity and integrity of blood-based evidence for accurate analysis.



FS754

SPECIFICATIONS

Performance

- Stabilizes blood-based prints and impressions
- Prevents smearing, distortion, and detail loss
- Prepares evidence for water-based enhancement methods
- Preserves clarity and contrast during processing
- Effective on porous and absorbent surfaces

Capacity

- 500 mL

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI