

Amido Black

Amido Black is a forensic reagent used to enhance latent fingerprints contaminated with blood by reacting with proteins in the residue. It produces a dark blue-black stain for clear visualization on a variety of surfaces.



SPECIFICATIONS

Performance

- Enhances blood-contaminated latent prints
- Reacts with proteins for strong, visible contrast
- Produces blue-black staining for clear detail
- Works on both porous and nonporous surfaces
- Not effective on normal (non-blood) latent prints

Variations

Model No.	Description	Capacity
FS762	Amido Black Powder	25 grams
FS763	Amido Black Methanol Base	500 mL
FS764	Amido Black Methanol Base	1L
FS765	Amido Black Methanol Rinse	500 mL
FS766	Amido Black Methanol Rinse	1L
FS767	Amido Black Aqueous	500 mL
FS768	Amido Black Aqueous	1L

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI