

# Roll-Ease Fingerprint Ink Rollers

These fingerprint ink rollers feature durable, unbreakable tapered handles designed for comfort and ease of use with any hand size. The stainless steel frame keeps the inked roller off surfaces when not in use, helping maintain cleanliness.



FS816-FS818

## SPECIFICATIONS

### Construction

- Tapered handle
- Frame
- Roller

### Options and Accessories

- Replacement rubber rollers (FS819-FS821)

### Performance

- Comfortable, tapered handles that fit any hand size
- The frame prevents the roller from contacting the tabletop when not in use

### Materials

- Rubber
- Stainless steel
- Plastic

## Variations

Model No.	Size
FS816	2"
FS817	3"
FS818	4"

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI