

Steel-Grip® Fingerprint Ink Rollers

These ink rollers feature smooth, solvent-resistant, medium-soft rubber for consistent ink application, paired with a durable stainless steel handle. The frame keeps the roller elevated when set down to prevent surface contact.



FS822 and FS823

SPECIFICATIONS

Construction

- Frame
- Roller

Dimensions

- Diameter: 1.1875"
- Handle: 5.25"

Performance

- The handle prevents the inked rubber roller from contacting the table surface

Materials

- Stainless steel
- Rubber

Options and Accessories

- Replacement fingerprint ink rubber rollers (FS819-FS821)

Variations

Model No.	Size
FS822	3"
FS823	4"

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI