


Scalpels - Retractable Blade Disposable

These spring-loaded, retractable safety scalpels are designed for safely collecting trace evidence from various surfaces. Available in multiple blade shapes, they offer controlled, precise scraping and can be easily cleaned between uses.



FS834-FS836

SPECIFICATIONS

Construction

- Spring-loaded
- 3 different blade shapes

Quantity

- 10 per package

Performance

- Great for scraping blood, paint, or other evidence from various surfaces
- Can be cleaned with 10% bleach solution or alcohol

Variations

Model No.	Description
FS834	Rounded #10
FS835	Rounded #15
FS836	Pointed #11

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI