


Latent Print Lifting Tape

This 3" x 180' tape is designed for collecting trace evidence such as hair, fibers, and glass fragments. Its water-soluble adhesive allows for easy recovery of collected evidence for further analysis.



FS916

SPECIFICATIONS

Performance

- For the collection of trace evidence such as hair, fiber, and glass shards
- Uses adhesive similar to fingerprint lifting
- Adhesive dissolves in water for easy evidence recovery
- Supports microscopic analysis after filtration

Dimensions

- 3" x 180'

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI